



Health And Hydrogen

HAH-1800

Flow Rate

H2 OUTPUT- 1200ml/min
O2 OUTPUT- 600 ml/min

FLOW RATE ADJUSTABLE
900/1200/1800/min

Advantages



One Click Run



Flow rate Adjustment



1-8 hours timing



Hydrogen only



Oxygen only



Oxy-hydrogen



Hydrogen water





Health And Hydrogen

HAH-1800



Handle

O₂(H₂) Outlet

H₂ Outlet

Water Control Panel

Smart Control Panel

Ventilation Holes

H₂ Window



Handle

Power Switch

Power plug

Water-level window

Thermovent

Product lable

Outfall




Water filling hole

Handle



Health And Hydrogen



One Click Run



Water Volume Monitoring



TDS Monitoring



1-8 Hours Time

*The diagrams in this manual are for reference only, please refer to the actual product

*If there is any technical change to this product in the future, such as the appearance and color of the product, it is subject to change without notice.

*This product is limited operation and must be used in accordance with the instruction manual

Hydrogen Accessories



NC-A silicon
Nasal Cannula



Y-S short
Y tube



Y-L Long
Y tube




NC-N nose
accessory



HR-A hydrogen
rod



Health And Hydrogen



H₂
O₂

HAH-1800

Product Parameters

Flow Rate

H₂ OUTPUT- 1200ml/min

O₂ OUTPUT- 600 ml/min

FLOW RATE ADJUSTABLE

900/1200/1800/min



Hydrogen purity

99.99%

Power input

AC220V/110V, 50Hz

Temperature

5°C~40°C

Pressure

700hPa-1060hPa

Relative Humidity

≤80%

Water quality requirement

Double distilled water or pure water

Product Size

300*193*285mm

Net Weight

6.5KG

Technology

SPE & PEM

Function

Hydrogen inhalation

Produce H₂ water with hydration rod

Advantage

High capacity / Easy Operation