







Health And Hydrogen

HAH-450

Flow Rate

H₂ OUTPUT- 300ml/min
O₂ OUTPUT- 150ml/min

Advantages

-  One Click Run
-  Flow rate Adjustment
-  1-8 hours timing
-  Hydrogen only
-  Oxygen only
-  Oxy-hydrogen
-  Hydrogen water





Health And Hydrogen

HAH-450



Handle

O₂(H₂) Outlet


H₂ Outlet

Water Control Panel

Smart Control Panel

Ventilation Holes

H₂ Window



Handle

Power Switch

Power plug

Water-level window

Thermovent

Product lable

Outfall




Water filling hole

Handle



Health And Hydrogen



One Click Run



Water Volume Monitoring



TDS Monitoring



1-8 Hours Time

*The diagrams in this manual are for reference only, please refer to the actual product

*If there is any technical change to this product in the future, such as the appearance and color of the product, it is subject to change without notice.

*This product is limited operation and must be used in accordance with the instruction manual

Hydrogen Accessories



NC-A silicon
Nasal Cannula



Y-S short
Y tube



Y-L Long
Y tube



NC-N nose
accessory



HR-A hydrogen
rod



Health And Hydrogen



HAH-450

Product Parameters

Flow Rate

H₂ OUTPUT- 300ml/min
O₂ OUTPUT- 150ml/min



Hydrogen purity

99.99%

Power input

AC220V/110V, 50Hz

Temperature

5°C~40°C

Pressure

700hPa-1060hPa

Relative Humidity

≤80%

Water quality requirement

Double distilled water or pure water

Product Size

300*193*285mm

Net Weight

6.5KG

Technology

SPE & PEM

Function

Hydrogen inhalation
Produce H₂ water with hydration rod

Advantage

High capacity / Easy Operation